



**R-PILOT™**  
**discover a new**  
**world of endo**

Your glide path file –  
if and when  
needed


# The first motor driven glide path instrument in reciprocating motion by VDW




- Glide path management can be a crucial step in root canal preparation
- With R-PILOT™, there finally is a safe and easy-to-use solution for cases where a glide path is needed
- R-PILOT™ is designed to be used in reciprocating motion with a designated drive system with original VDW RECIPROC® settings

## Designed for reciprocation – optimized for easy usage


- One size
- Length: 21 mm, 25 mm, 31 mm
- Material: M-wire
- Constant taper: 4 %
- Single use



Cutting efficiency  
with S-shape cross section



# Easy and efficient glide path management with R-PILOT™



## Safe

- Less risk of instrument fractures\*
- Reduction of cyclic fatigue\*
- Single patient use instrument



## Predictable shape

- One file glide path management
- Better respect of the canal anatomy\*\*



## Time saving

- Sterile instruments and less files needed compared to manual files
- No need to change motor settings if you are already working with reciprocating motion

## R-PILOT™ Product range

R-PILOT™ NiTi Glidepath Files	Blister à 6 instruments	One size 12.5
1 x Blister of R-PILOT™ Glidepath Files	Length 21 mm	V04 0217 021 012
1 x Blister of R-PILOT™ Glidepath Files	Length 25 mm	V04 0217 025 012
1 x Blister of R-PILOT™ Glidepath Files	Length 31 mm	V04 0217 031 012

\*compared to rotary multiple use instruments  
\*\*versus VDW SST files





For further information please contact your preferred distributor or visit [www.vdw-dental.com](http://www.vdw-dental.com)

VDW GmbH  
Bayerwaldstr. 15  
81737 Munich  
Germany  
Tel +49 (0) 89 62734-0  
Fax +49 (0) 89 62734-304  
[info@vdw-dental.com](mailto:info@vdw-dental.com)

